

Die Line # 222 (Thick Canister Wrap around label)

4.25"

.75"

**Special Instructions for Cleaning and Decontamination Against HIV-1 and HBV of Surfaces or Objects Soiled With Blood or Body Fluids:**

**Personal Protection:** When using this disinfecting towelette, wear appropriate barrier protection such as disposable protective gloves, protective gown, face mask and eye covering.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfecting towelette. A Z<sup>3</sup> Wipes towelette can be used for this purpose.

**Contact Time:** Thoroughly wet surface with a new Z<sup>3</sup> Wipes towelette. Repeated use of the product may be required to ensure that the surface remains visibly wet for 2 minutes. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

**Disposal Of Infectious Materials:** Cleaning materials used that may contain blood or other body fluids should be autoclaved and disposed of in accordance with Federal, State and local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

**PRECAUTIONARY STATEMENTS**  
**Hazards to Humans and Domestic Animals**  
**CAUTION**

Harmful if absorbed through the skin. Causes moderate eye irritation. Avoid contact with eyes, skin, or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse.

**FIRST AID**

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for further advice.

**If on skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

For emergency information on product, call Customer Service at 1-800-841-1428, Monday through Friday, between 6 a.m. and 4 p.m. Pacific Time. After 4 p.m., call Infotrac 24-Hour Emergency Service at 1-800-535-5053.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

**STORAGE AND DISPOSAL**

**Towelette** - Do not reuse towelette. Dispose of used towelette in trash.  
**Dispenser** - Do not reuse empty container. Wrap container and discard in trash.

For product or technical information, contact Benco Dental at 1-800-GO-BENCO or visit our website at www.benco.com.

Distributed By:  
 Benco Dental Supply Co.  
 11 Bear Creek Blvd., Wilkes-Barre, PA 18702

70-1190-2  
 84604-22A

Reorder No:



[3738-967]

Z<sup>3</sup> Wipes Disinfecting Towelettes are non-woven disposable towelettes intended for use on hard, non-porous surfaces in health care settings such as dental offices, hospitals, emergency medical settings, laboratories and clinics. Z<sup>3</sup> Wipes are effective against the following microorganisms when used as directed:

- *Mycobacterium tuberculosis* var: *bovis* (BCG)
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella choleraesuis*
- *Trichophyton mentagrophytes*
- Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Vancomycin Resistant *Enterococcus faecalis* (VRE)
- Hepatitis B Virus (HBV)\*
- Herpes Simplex Virus Types 1 and 2\*
- Human Immunodeficiency Virus (HIV-1)\*



- Tuberculocidal
- Bactericidal
- Virucidal\*
- Fungicidal

**ACTIVE INGREDIENTS:**

Diisobutylphenoxethoxyethyl dimethyl benzyl ammonium chloride.....	0.28%
Isopropanol.....	17.20%
INERT INGREDIENTS†.....	82.52%
<b>TOTAL</b> .....	<b>100.00%</b>

† Does not include the weight of the dry towelette.

**KEEP OUT OF REACH OF CHILDREN**

**CAUTION** See other panel for additional precautionary statements.

EPA Reg. No. 46781-8-35659  
 EPA Est. No. 56952-WI-001

Made in the U.S.A.



NET CONTENTS: 1 lb. 13 oz. • 160 Pre-Saturated Towelettes (6" x 6.76")

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Z<sup>3</sup> Wipes are non-woven disposable towelettes intended for use on hard, non-porous surfaces in health care settings such as dental offices, hospitals, emergency medical settings, laboratories, clinics, and other critical care areas where control of cross contamination between treated surfaces is required.

Z<sup>3</sup> Wipes will clean and disinfect, when used as directed, such items as operating room tables and lights, laboratory equipment and surfaces, dental unit instrument trays, operatory light switches, exterior surfaces of amalgamators and dental curing lights, exterior surfaces of anesthesia machines and respiratory therapy equipment, dental countertops, dental chairs, loupes, exterior surfaces of endodontic equipment, such as apex locators, pulp testers and motors, stools, floors, walls, handrails, door knobs, tables, chairs, faucets, bathroom fixtures, exterior toilet surfaces, garbage cans, computers, keyboards, cabinets, lights and telephones, which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.

**DO NOT USE AS A BABY WIPE.**

**Not For Use On Skin**

**Dispenser Directions:**

To start feed, remove cover and discard seal. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at a 90° angle. Remaining towelettes

feed automatically, ready for next use. When not in use, keep center cap closed to prevent solution loss.

**Cleaning Instructions:**

Use one Z<sup>3</sup> Wipes towelette to completely preclean surface of all gross debris.

**For use as a virucide\*, including HIV-1 and HBV, Methicillin Resistant *Staphylococcus aureus* (MRSA) and Vancomycin Resistant *Enterococcus faecalis* (VRE):**

Use a second Z<sup>3</sup> Wipes towelette to thoroughly wet the surface. Repeated use of the product may be required to ensure that the surface remains visibly wet for 2 minutes at room temperature (69 °F/20 °C).

**For use as a disinfectant:**

Use a second Z<sup>3</sup> Wipes towelette to thoroughly wet the surface. Repeated use of the product may be required to ensure that the surface remains visibly wet for 3 minutes at room temperature (69 °F/20 °C) for Tuberculocidal activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella choleraesuis*.

Z<sup>3</sup> Wipes towelettes effectively kills Human Immunodeficiency Virus Type 1 (HIV) and Hepatitis B Virus (HBV) on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HIV and HBV.

No Print Area • Glue Adhesive

16.75"

- Critical Text Box
- - - - - Bleeds
- - - - - Adhesive area
- \_\_\_\_\_ Die Cut area