

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a Poison Control Center or doctor for treatment advice.

IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a Poison Control Center or doctor for treatment advice.

IF SWALLOWED: Call a Poison Control Center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told by a Poison Control Center or doctor. Do not give anything by mouth to an unconscious person.

IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration preferably mouth-to-mouth if possible. Call a Poison Control Center or doctor for further treatment advice.

In case of an emergency for additional information call toll free 800-535-5053. Have the product container or label with you when calling a Poison Control Center or doctor, or going for treatment.

NOTE TO PHYSICIAN
Probable mucosal damage may contraindicate the use of gastric lavage.

MANUFACTURED BY:
CERTOL
INTERNATIONAL, LLC
6120 E. 58TH AVE.
COMMERCE CITY, CO 80022 USA
1 800 THE EDGE
(1-800-843-3343)

ProSpray C-60™

Concentrated Surface Disinfectant and Holding Solution

Lemon Scent

Tuberculocidal, Virucidal, Bactericidal, Fungicidal

Active Ingredients

o-phenylphenol	9.0%
o-benzyl-p-chlorophenol	1.0%
Inert Ingredients	90.0%
Total	100.0%

E. P.A. Reg No. 46851-1
E. P.A. Est. No. 46851-CO-003

**Contents: 1 Quart
(32 fluid oz. - 946 ml.)
Reorder No. PSC60/32**

Lot Number is printed on container

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV ON SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS

PERSONAL PROTECTION: When handling items soiled with blood or body fluids or when conducting procedures which might lead to contact with blood or body fluids, always wear items such as protective gloves, gowns, masks and eye coverings.

CLEANING PROCEDURES: Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: ProSpray C-60 when diluted 1:32 has been shown to be effective in destroying the HIV-1 (AIDS virus) in 60 seconds at 20°C. This contact time will not control other common types of viruses and bacteria.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

PRECAUTIONARY STATEMENT HAZARD TO HUMANS AND DOMESTIC ANIMALS

DANGER: Corrosive. Causes irreversible eye damage and skin burns. Harmful if inhaled, swallowed or absorbed through skin. Do not get in eyes, skin or clothing. Avoid breathing spray mist. Wear goggles or face shields (safety glasses), protective clothing (long-sleeve shirt and long pants, socks, plus shoes and chemical resistant gloves such as water proof gloves).

USER SAFETY REQUIREMENTS: Follow manufacturers' instructions for cleaning/maintaining PPE. If no such instructions for washables, use detergent and hot water. Keep and wash PPE separately from other laundry.

USER SAFETY RECOMMENDATIONS: Users must wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Users must remove clothing immediately if pesticide gets inside. Then wash thoroughly and put on clean clothing. Users must remove PPE immediately after handling this product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and change into clean clothing.

DIRECTIONS FOR USE IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

ProSpray C-60 is designed for use on hard non-porous surfaces and/or in holding solutions in medical, dental, laboratory and veterinary facilities. Thoroughly clean surface of items to be disinfected and rough dry. ProSpray C-60 is a concentrated solution and must be diluted 1:32 with water prior to instrument soak and/or decontamination of dental and medical instruments or disinfection of hard environmental surfaces. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Soak/Decontamination – immerse items in diluted ProSpray C-60 solution for instrument soak and/or decontamination of medical and dental instruments. Rinse articles thoroughly following soak/decontamination. Organic materials or water introduced into the diluted product may effect products effectiveness. Discard the used solution daily or sooner if it becomes diluted or soiled.

Ultrasonic Cleaning – Place instruments into diluted ProSpray C-60 and activate ultrasonic unit. When instruments are clean, remove and rinse thoroughly. Change solution daily or sooner if it becomes diluted or soiled.

Cleaning – Apply ProSpray C-60 directly to surface using mop or spray bottle. Wipe surface to remove debris and solution using a clean paper or cloth towel and discard.

Hard Surface Disinfection – diluted ProSpray C-60 solution should be mopped, sponged, sprayed or wiped onto hard environmental surfaces such as walls, floors and countertops in a manner to thoroughly wet the surface. A properly prepared solution of ProSpray C-60 intended for use as a hard surface disinfectant, has a shelf life of 60 days (use life) when stored in a closed container such as a spray bottle.

A contact time of 10 minutes at 20°C (68°F) with a 1:32 dilution of ProSpray C-60 has been shown to be effective on bacteria including: Escherichia Coli 0113, Staphylococcus aureus, Salmonella typhosa, Pseudomonas aeruginosa, Mycobacterium tuberculosis var.bovis; fungi including: Trichophyton mentagrophytes; viruses including Influenza virus type A1 (New Jersey) and Influenza B (Maryland), Herpes virus types 1 and 2, and Vaccinia virus.*

ProSpray C-60 when diluted 1:32 has been shown to be effective in destroying the HIV-1 (AIDS virus) in 60 seconds at 20°C (68°F). Although efficacy in 60 second contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use 10 minute contact time for disinfection against all the organisms claimed. Effective in waters up to 400 ppm hardness calculated as calcium carbonate.

For use on metals, plastics, rubber and glass. Glass may etch over prolonged or repeated exposure to

concentrated solutions. Do not store lensed articles in concentrated ProSpray C-60.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Keep ProSpray C-60 covered in storage or in use.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL: Do not reuse empty container. Triple rinse container, puncture and dispose of in a sanitary landfill, or incinerate.



+H816PSC60-32D

TO PREPARE 1:32 USE-DILUTION

	To Make: 1 Pint	To Make: 1 Quart	To Make: 1/2 Gallon	To Make: 1 Gallon
Water	15 1/2 fl. oz.	31 fl. oz.	62 fl. oz.	124 fl. oz.
ProSpray C-60	1/2 fl. oz. (15 ml.)	1 fl. oz. (30 ml.)	2 fl. oz. (60 ml.)	4 fl. oz. (120 ml.)

